




MIS Pedicle Screw System



Bringing Simplicity and Elegance to Minimally Invasive Surgery

MIS Pedicle Screw System



Removable Tension Ring


Keeps tab secure during placement and manipulation of screw, easily removed for ease in removing elongated tabs

Integrated Tower

Break-off tabs allow for less soft tissue disruption


Polyaxial Screw

70° conical angulation provides flexibility



Superior MIS Rods

Distal tip of rod aids in penetrating tissue, proximal face allows for full control of rod trajectory



Integrated Reduction Threads

Aid in securing rod in complex anatomies




Self-Tapping Threads

More secure bone engagement



Dual Lead Threads

Provides quick and efficient screw placement



Cannulated Screw

Available in both solid and cannulated screws

