




Pedicle Screw System




Multifaceted System that Adapts to Meet Multiple Pathologies

Pedicle Screw System

Polyaxial Screw

70° conical angulation provides flexibility



Full Complement of Implants

Includes polyaxial, reduction, and mono screws, with cortical and bi-cortical options. Hooks, offset connectors, multiple rod to rod connectors, and single joint and dual joint crosslinks

Curved and Straight Rods

Rods had etched line for better visualization and hex drives to aid in manipulation



Dual Lead Threads

Provides quick and efficient screw placement

Cannulated Screw

Available in both solid and cannulated screws



Self-Tapping Threads

More secure bone engagement



Hexalobe Drives

Both the bone screw and set screw have large hexalobe drives (T25 and T40) for superior control and torque

